

## #29-092 – Edmund Optics® Lab Storage Kit



Edmund Optics® Lab Storage Kit provide an ideal solution for protecting and organizing components in an optics lab environment. Each Lab Storage Kit includes 24 component storage drawers, 12 optics storage trays, drawer labels, and a marker. These convenient storage cases can be placed on a lab bench, easily wall mounted, and will help alleviate clutter and prevent loss and damage to precision components

### Specifications

<b>Cabinet Material:</b>	Polystyrene	<b>Includes Dividers:</b>	<b>Yes</b>
<b>Cabinet Color:</b>	Black	<b>Mounting Location:</b>	Wall
<b>Cabinet Style:</b>	Drawer	<b>Number of Dividers:</b>	12
<b>Drawer Color:</b>	Clear	<b>Number of Drawers:</b>	24
<b>Drawer Inside Depth:</b>	4-9/16 in	<b>Overall Depth:</b>	6-1/2 in
<b>Drawer Inside Height HxWxD</b>	6" x 4-9/16" x 2-3/16"	<b>Overall Height:</b>	15-3/4 in
<b>Drawer Inside Width:</b>	6 in	<b>Overall Width:</b>	20 in
<b>Drawer Material:</b>	Plastic	<b>Stackable:</b>	Yes